



MADCAM1 Antibody, Biotin conjugated

Product Code	CSB-PA618776HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13477
Immunogen	Recombinant Human Mucosal addressin cell adhesion molecule 1 protein (19-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mucosal addressin cell adhesion molecule 1 (MAdCAM-1) (hMAdCAM-1), MADCAM1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MADCAM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.