



GTF2H2 Antibody

Product Code	CSB-PA618808LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13888
Immunogen	Recombinant Human General transcription factor IIH subunit 2 protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	General transcription factor IIH subunit 2 (Basic transcription factor 2 44 kDa subunit) (BTF2 p44) (General transcription factor IIH polypeptide 2) (TFIIH basal transcription factor complex p44 subunit), GTF2H2, BTF2P44
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GTF2H2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.