




LRRC32 Antibody

Product Code	CSB-PA618886LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14392
Immunogen	Recombinant Human Transforming growth factor beta activator LRRC32 protein (301-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta activator LRRC32 (Garpin) (Glycoprotein A repetitions predominant) (GARP) (Leucine-rich repeat-containing protein 32), LRRC32, D11S833E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LRRC32

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA618886LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.