



LRRC32 Antibody, HRP conjugated

Product Code	CSB-PA618886LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14392
Immunogen	Recombinant Human Transforming growth factor beta activator LRRC32 protein (301-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta activator LRRC32 (Garpin) (Glycoprotein A repetitions predominant) (GARP) (Leucine-rich repeat-containing protein 32), LRRC32, D11S833E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LRRC32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.