



HERC3 Antibody, Biotin conjugated

Product Code	CSB-PA618972LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15034
Immunogen	Recombinant Human Probable E3 ubiquitin-protein ligase HERC3 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable E3 ubiquitin-protein ligase HERC3 (EC 2.3.2.26) (HECT domain and RCC1-like domain-containing protein 3) (HECT-type E3 ubiquitin transferase HERC3), HERC3, KIAA0032
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	HERC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.