



RHOH Antibody, HRP conjugated

Product Code	CSB-PA618992LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15669
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoH protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoH (GTP-binding protein TTF) (Translocation three four protein), RHOH, ARHH TTF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.