



RHOH Antibody, HRP conjugated

Product Code	CSB-PA618992LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15669
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoH protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoH (GTP-binding protein TTF) (Translocation three four protein), RHOH, ARHH TTF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOH