



PCK2 Antibody, HRP conjugated

Product Code	CSB-PA619092LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16822
Immunogen	Recombinant Human Phosphoenolpyruvate carboxykinase [GTP], mitochondrial protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoenolpyruvate carboxykinase [GTP], mitochondrial (PEPCK-M) (EC 4.1.1.32), PCK2, PEPCK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PCK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.