



PSMD2 Antibody, HRP conjugated

Product Code	CSB-PA619755LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13200
Immunogen	Recombinant Human 26S proteasome non-ATPase regulatory subunit 2 protein (401-700AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome non-ATPase regulatory subunit 2 (26S proteasome regulatory subunit RPN1) (26S proteasome regulatory subunit S2) (26S proteasome subunit p97) (Protein 55.11) (Tumor necrosis factor type 1 receptor-associated protein 2), PSMD2, TRAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.