



SMAD4 Antibody, FITC conjugated

Product Code	CSB-PA619768LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13485
Immunogen	Recombinant Human Mothers against decapentaplegic homolog 4 protein (27-195AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mothers against decapentaplegic homolog 4 (MAD homolog 4) (Mothers against DPP homolog 4) (Deletion target in pancreatic carcinoma 4) (SMAD family member 4) (SMAD 4) (Smad4) (hSMAD4), SMAD4, DPC4 MADH4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SMAD4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.