




KRT72 Antibody

Product Code	CSB-PA619911LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14CN4
Immunogen	Recombinant Human Keratin, type II cytoskeletal 72 protein (301-511AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 72 (Cytokeratin-72) (CK-72) (Keratin-72) (K72) (Type II inner root sheath-specific keratin-K6irs2) (Type-II keratin Kb35), KRT72, K6IRS2 KB35 KRT6 KRT6IRS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT72

Image



Western blot
 All lanes: KRT72 antibody at 2µg/ml + MCF-7 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 56, 52 kDa
 Observed band size: 56, 52 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.