



KRT72 Antibody, Biotin conjugated

Product Code	CSB-PA619911LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14CN4
Immunogen	Recombinant Human Keratin, type II cytoskeletal 72 protein (301-511AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 72 (Cytokeratin-72) (CK-72) (Keratin-72) (K72) (Type II inner root sheath-specific keratin-K6irs2) (Type-II keratin Kb35), KRT72, K6IRS2 KB35 KRT6 KRT6IRS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT72
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.