



MAPRE1 Antibody, Biotin conjugated

Product Code	CSB-PA620890LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15691
Immunogen	Recombinant Human Microtubule-associated protein RP/EB family member 1 protein (2-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein RP/EB family member 1 (APC-binding protein EB1) (End-binding protein 1) (EB1), MAPRE1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MAPRE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.