



DNASE1L3 Antibody, Biotin conjugated

Product Code	CSB-PA621686ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13609
Immunogen	Recombinant Human Deoxyribonuclease gamma protein (21-305AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribonuclease gamma (DNase gamma) (EC 3.1.21.-) (DNase I homolog protein DHP2) (Deoxyribonuclease I-like 3) (DNase I-like 3) (Liver and spleen DNase) (LS-DNase) (LSD), DNASE1L3, DHP2 DNAS1L3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DNASE1L3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.