



TUBA3C Antibody, HRP conjugated

Product Code	CSB-PA621695LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13748
Immunogen	Recombinant Human Tubulin alpha-3C/D chain protein (161-418AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin alpha-3C/D chain (Alpha-tubulin 2) (Alpha-tubulin 3C/D) (Tubulin alpha-2 chain) [Cleaved into: Detyrosinated tubulin alpha-3C/D chain], TUBA3C; TUBA3D, TUBA2;
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TUBA3C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.