



INSL4 Antibody, HRP conjugated

Product Code	CSB-PA621784LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14641
Immunogen	Recombinant Human Early placenta insulin-like peptide protein (30-139AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Early placenta insulin-like peptide (EPIL) (Insulin-like peptide 4) (Placentin) [Cleaved into: Early placenta insulin-like peptide B chain; Early placenta insulin-like peptide A chain], INSL4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INSL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.