



RTL8B Antibody, HRP conjugated

Product Code	CSB-PA622040LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17RB0
Immunogen	Recombinant Human Retrotransposon Gag-like protein 8B protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retrotransposon Gag-like protein 8B (Mammalian retrotransposon derived protein 8B), RTL8B, CXX1C FAM127C MAR8B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTL8B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.