



RTL8B Antibody, FITC conjugated

Product Code	CSB-PA622040LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17RB0
Immunogen	Recombinant Human Retrotransposon Gag-like protein 8B protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retrotransposon Gag-like protein 8B (Mammalian retrotransposon derived protein 8B), RTL8B, CXX1C FAM127C MAR8B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTL8B