



SLC10A2 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA622519LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q12908 |
| Immunogen | Recombinant Human Ileal sodium/bile acid cotransporter protein (306-348AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ileal sodium/bile acid cotransporter (Apical sodium-dependent bile acid transporter) (ASBT) (Ileal Na(+)/bile acid cotransporter) (Ileal sodium-dependent bile acid transporter) (IBAT) (ISBT) (Na(+)-dependent ileal bile acid transporter) (Sodium/taurocholate cotransporting polypeptide, ileal) (Solute carrier family 10 member 2), SLC10A2, ASBT ISBT NTCP2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | SLC10A2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |