



SLC10A2 Antibody, FITC conjugated

Product Code	CSB-PA622519LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12908
Immunogen	Recombinant Human Ileal sodium/bile acid cotransporter protein (306-348AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ileal sodium/bile acid cotransporter (Apical sodium-dependent bile acid transporter) (ASBT) (Ileal Na(+)/bile acid cotransporter) (Ileal sodium-dependent bile acid transporter) (IBAT) (ISBT) (Na(+)-dependent ileal bile acid transporter) (Sodium/taurocholate cotransporting polypeptide, ileal) (Solute carrier family 10 member 2), SLC10A2, ASBT ISBT NTCP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC10A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.