





LAPTM5 Antibody

Product Code	CSB-PA622667LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13571
Immunogen	Recombinant Human Lysosomal-associated transmembrane protein 5 protein (205-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosomal-associated transmembrane protein 5 (Lysosomal-associated multitransmembrane protein 5) (Retinoic acid-inducible E3 protein), LAPTM5, KIAA0085
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	LAPTM5

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA622667LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA622667LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.