



KCNJ11 Antibody, Biotin conjugated

Product Code	CSB-PA622781LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14654
Immunogen	Recombinant Human ATP-sensitive inward rectifier potassium channel 11 protein (172-390AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-sensitive inward rectifier potassium channel 11 (IKATP) (Inward rectifier K(+) channel Kir6.2) (Potassium channel, inwardly rectifying subfamily J member 11), KCNJ11
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNJ11