




GABRA6 Antibody

Product Code	CSB-PA623078LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16445
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor subunit alpha-6 protein (20-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid receptor subunit alpha-6 (GABA(A) receptor subunit alpha-6), GABRA6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GABRA6

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA623078LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.