



GABRA6 Antibody, Biotin conjugated

Product Code	CSB-PA623078LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16445
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor subunit alpha-6 protein (20-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid receptor subunit alpha-6 (GABA(A) receptor subunit alpha-6), GABRA6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GABRA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.