



FAM187B Antibody, HRP conjugated

Product Code	CSB-PA623168LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17R55
Immunogen	Recombinant Human Protein FAM187B protein (18-335AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM187B (Transmembrane protein 162), FAM187B, TMEM162
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM187B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.