





MAB21L1 Antibody

Product Code	CSB-PA623800LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13394
Immunogen	Recombinant Human Putative nucleotidyltransferase MAB21L1 protein (1-359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative nucleotidyltransferase MAB21L1 (EC 2.7.7.-) (Protein mab-21-like 1), MAB21L1, CAGR1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAB21L1

Image



Western blot
 All lanes: MAB21L1 antibody at 6µg/ml
 Lane 1: A375 whole cell lysate
 Lane 2: CEM whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 41 kDa
 Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA623800LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.