



MAB21L1 Antibody, Biotin conjugated

Product Code	CSB-PA623800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13394
Immunogen	Recombinant Human Putative nucleotidyltransferase MAB21L1 protein (1-359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative nucleotidyltransferase MAB21L1 (EC 2.7.7.-) (Protein mab-21-like 1), MAB21L1, CAGR1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAB21L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.