



SMCO1 Antibody, HRP conjugated

Product Code	CSB-PA623935LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q147U7
Immunogen	Recombinant Human Single-pass membrane and coiled-coil domain-containing protein 1 protein (82-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-pass membrane and coiled-coil domain-containing protein 1 (Single-pass membrane protein with coiled-coil domains 1), SMCO1, C3orf43
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SMCO1