



# KCTD19 Antibody, FITC conjugated

Product Code	CSB-PA624192LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17RG1
Immunogen	Recombinant Human BTB/POZ domain-containing protein KCTD19 protein (490-751AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing protein KCTD19, KCTD19
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCTD19