



NGF Antibody

Product Code	CSB-PA643662LA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29074
Immunogen	Recombinant Pig Beta-nerve growth factor protein (110-229AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-nerve growth factor (Beta-NGF) (Fragment), NGF, NGFB
Immunogen Species	Sus scrofa (Pig)
Research Area	Others
Target Names	NGF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.