





HKDC1 Antibody

Product Code	CSB-PA647983LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2TB90
Immunogen	Recombinant Human Putative hexokinase HKDC1 protein (129-427AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative hexokinase HKDC1 (EC 2.7.1.1) (Hexokinase domain-containing protein 1), HKDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HKDC1

Image



Western blot
All lanes: HKDC1 antibody at 4µg/ml + U87 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 103, 91, 83 kDa
Observed band size: 103 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA647983LA01HU at dilution of 1:100