



HKDC1 Antibody, Biotin conjugated

Product Code	CSB-PA647983LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2TB90
Immunogen	Recombinant Human Putative hexokinase HKDC1 protein (129-427AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative hexokinase HKDC1 (EC 2.7.1.1) (Hexokinase domain-containing protein 1), HKDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HKDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.