



ZXDC Antibody, Biotin conjugated

Product Code	CSB-PA648494LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2QGD7
Immunogen	Recombinant Human Zinc finger protein ZXDC protein (139-416AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein ZXDC (ZXD-like zinc finger protein), ZXDC, ZXDL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ZXDC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.