



KLHL38 Antibody, FITC conjugated

Product Code	CSB-PA652063LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2WGJ6
Immunogen	Recombinant Human Kelch-like protein 38 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kelch-like protein 38, KLHL38, C8orfK36
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KLHL38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.