



PRICKLE4 Antibody, Biotin conjugated

Product Code	CSB-PA652689LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2TBC4
Immunogen	Recombinant Human Prickle-like protein 4 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prickle-like protein 4 (Overexpressed breast tumor protein), PRICKLE4, C6orf49 OEBT
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRICKLE4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.