



PGS1 Antibody, Biotin conjugated

Product Code	CSB-PA653479LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q32NB8
Immunogen	Recombinant Human CDP-diacylglycerol--glycerol-3-phosphate 3-phosphatidyltransferase, mitochondrial protein (110-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CDP-diacylglycerol--glycerol-3-phosphate 3-phosphatidyltransferase, mitochondrial (EC 2.7.8.5) (Phosphatidylglycerophosphate synthase 1) (PGP synthase 1), PGS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PGS1