



KRTAP19-5 Antibody

Product Code	CSB-PA659465LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3LI72
Immunogen	Recombinant Human Keratin-associated protein 19-5 protein (1-72AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 19-5, KRTAP19-5, KAP19.5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KRTAP19-5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.