



TSGA10IP Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA660950LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q3SY00 |
| Immunogen | Recombinant Human Testis-specific protein 10-interacting protein (6-302AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Testis-specific protein 10-interacting protein (Tsga10-interacting protein), TSGA10IP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | TSGA10IP |