



TSGA10IP Antibody, Biotin conjugated

Product Code	CSB-PA660950LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY00
Immunogen	Recombinant Human Testis-specific protein 10-interacting protein (6-302AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testis-specific protein 10-interacting protein (Tsga10-interacting protein), TSGA10IP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TSGA10IP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.