





KRT71 Antibody

Product Code	CSB-PA660951LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY84
Immunogen	Recombinant Human Keratin, type II cytoskeletal 71 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 71 (Cytokeratin-71) (CK-71) (Keratin-71) (K71) (Type II inner root sheath-specific keratin-K6irs1) (Keratin 6 irs) (hK6irs) (hK6irs1) (Type-II keratin Kb34), KRT71, K6IRS1 KB34 KRT6IRS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT71

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA660951LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA660951LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.