



KRT71 Antibody, HRP conjugated

Product Code	CSB-PA660951LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY84
Immunogen	Recombinant Human Keratin, type II cytoskeletal 71 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type II cytoskeletal 71 (Cytokeratin-71) (CK-71) (Keratin-71) (K71) (Type II inner root sheath-specific keratin-K6irs1) (Keratin 6 irs) (hK6irs) (hK6irs1) (Type-II keratin Kb34), KRT71, K6IRS1 KB34 KRT6IRS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT71
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.