





FAM161A Antibody

Product Code	CSB-PA661553LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3B820
Immunogen	Recombinant Human Protein FAM161A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM161A, FAM161A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM161A


Image



Western blot
 All lanes: FAM161A antibody at 8µg/ml
 Lane 1: Human placenta tissue
 Lane 2: Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 77, 65, 84 kDa
 Observed band size: 77 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA661553LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA661553LA01HU at dilution of 1:100