




# FAM161A Antibody

<b>Product Code</b>	CSB-PA661553LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q3B820
<b>Immunogen</b>	Recombinant Human Protein FAM161A protein (1-300AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Protein FAM161A, FAM161A
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	FAM161A

## Image



### Western blot

All lanes: FAM161A antibody at 8µg/ml

Lane 1: Human placenta tissue


Lane 2: Mouse liver tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 77, 65, 84 kDa

Observed band size: 77 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA661553LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA661553LA01HU at dilution of 1:100

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.