



FAM161A Antibody, Biotin conjugated

Product Code	CSB-PA661553LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3B820
Immunogen	Recombinant Human Protein FAM161A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM161A, FAM161A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM161A