



GRAMD1B Antibody, Biotin conjugated

Product Code	CSB-PA662055LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3KR37
Immunogen	Recombinant Human GRAM domain-containing protein 1B protein (235-504AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GRAM domain-containing protein 1B, GRAMD1B, KIAA1201
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GRAMD1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.