



SLC25A52 Antibody, Biotin conjugated

Product Code	CSB-PA662246LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SY17
Immunogen	Recombinant Human Solute carrier family 25 member 52 protein (1-33AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 25 member 52 (Mitochondrial carrier triple repeat protein 2), SLC25A52, MCART2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A52
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.