



# MOGAT2 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA663570LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q3SYC2
<b>Immunogen</b>	Recombinant Human 2-acylglycerol O-acyltransferase 2 protein (124-284AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	2-acylglycerol O-acyltransferase 2 (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 2) (MGAT2) (hMGAT2) (Diacylglycerol O-acyltransferase candidate 5) (hDC5) (Diacylglycerol acyltransferase 2-like protein 5) (Monoacylglycerol O-acyltransferase 2), MOGAT2, DC5 DGAT2L5
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	MOGAT2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.