



MOGAT2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA663570LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q3SYC2 |
| Immunogen | Recombinant Human 2-acylglycerol O-acyltransferase 2 protein (124-284AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 2-acylglycerol O-acyltransferase 2 (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 2) (MGAT2) (hMGAT2) (Diacylglycerol O-acyltransferase candidate 5) (hDC5) (Diacylglycerol acyltransferase 2-like protein 5) (Monoacylglycerol O-acyltransferase 2), MOGAT2, DC5 DGAT2L5 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | MOGAT2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |