



MOGAT2 Antibody, Biotin conjugated

Product Code	CSB-PA663570LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3SYC2
Immunogen	Recombinant Human 2-acylglycerol O-acyltransferase 2 protein (124-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	2-acylglycerol O-acyltransferase 2 (EC 2.3.1.22) (Acyl-CoA:monoacylglycerol acyltransferase 2) (MGAT2) (hMGAT2) (Diacylglycerol O-acyltransferase candidate 5) (hDC5) (Diacylglycerol acyltransferase 2-like protein 5) (Monoacylglycerol O-acyltransferase 2), MOGAT2, DC5 DGAT2L5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MOGAT2