



# SLC10A6 Antibody, HRP conjugated

Product Code	CSB-PA665974LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3KNW5
Immunogen	Recombinant Human Solute carrier family 10 member 6 protein (312-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 10 member 6 (Sodium-dependent organic anion transporter), SLC10A6, SOAT
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC10A6